User Manual

MY01 Mobile Application

MY01 Inc. 400 Boul de Maisonneuve Ouest, Suite 700, Montréal, QC, Canada, H3A 1L4 +1 (855) 799-6901

MedEnvoy Prinses Margrietplantsoen 33 Suite 123, 2595 AM The Hague The Netherlands

CEpartner4U BV Esdoornlaan 13 3951 DB Maarn The Netherlands

Introduction	3
Getting Started	4
Recommended smart device settings	4
Bluetooth	4
Date & Time	5
Signing up	6
Signing in	6
Application Interface	6
Settings	8
Instructions	9
Connecting to the MY01 device	9
Active Device information	11
Remote Data Display	11
Disconnecting to the MY01 device	13
List of inactive devices	14
Logging out of the MY01 Mobile Application	15
Data retrieval	16
Help Center/ FAQs	17
Chat	18
Technical Issues	18

Symbol Explanation	Symbol
Caution : Federal law (USA) restricts this device to sale by or on the order of a physician	R_x ONLY
Manufacturer	
Authorized representative in the European Community	EC REP
Medical Device	MD
Unique Device Identifier (UDI)	UDI
Importer	

Introduction

Intended use

The MY01 Mobile Application is an optional application intended for storing and displaying identical pressure values from the MY01 Continuous Compartmental Pressure Monitor device. The data is for informational purposes only and is not intended to be used for diagnosis of any nature or active patient monitoring.

Patient target group

The device is intended for use in patients at risk or with suspected compartment syndrome due to a wide variety of clinical presentations including trauma.

Intended Users

Trained trauma surgeons/Physicians.

Application download

Getting Started

Recommended smart device settings

Ensure your smart device settings correspond with the following parameters:

<u>**Compatibility**</u>: For a list of smart devices and operating systems that work with the MY01 Mobile Application check <u>help.my01.io</u>

- Before downloading and installing the app, make sure your mobile device meets or exceeds the minimum performance specifications. Additional information on the most up to date device requirements can be found at <u>help.my01.io</u>
- The MY01 device uses Bluetooth Low Energy (BLE) technology specification 4.2 Ensure that the mobile phone uses BLE-compatible wireless technology.
- Recommended: Access to your mobile device should be password protected.

Bluetooth

Enabling Bluetooth on your iOS device

- 1. Tap on **Settings** on your iPhone.
- 2. Select General.
- 3. Tap on Bluetooth.
- 4. Toggle Bluetooth to **On.**

Enabling Bluetooth on your Android device:

- 1. Tap on **Settings** on your Android device.
- 2. Look for **Bluetooth** or the Bluetooth symbol in your settings and tap it.
- There should be an option to enable. Please tap or swipe on it so that is in the on position.

Date & Time

Update Date & Time on Your iOS Device

- 1. From your home screen, navigate to Settings.
- 2. Tap General.
- 3. Tap Date & Time.
- 4. Make sure that the option Set Automatically is turned on.
- 5. If this option is turned off, check that the correct Date, Time and Time Zone are selected.

Update Date & Time on Your Android Device

- 1. Tap Settings to open the Settings menu.
- 2. Tap Date & Time.
- 3. Tap Automatic.
- 4. If this option is turned off, check that the correct Date, Time and Time Zone are selected.

Signing up

To connect the MY01 device through Bluetooth to the MY01 Mobile Application for the first time, users must follow these steps:

- 1. Users are invited by their institution, and should follow the steps in the received email. Choose a secure password and confirm your account.
- Follow the link to install the MY01 Mobile Application from the Google Play or Apple AppStore.

Signing in

- 1. Open the MY01 Mobile Application.
- 2. Login using the same credentials as entered in the web portal.
- 3. If you've forgotten your password, click on reset password and follow the on-screen instructions.

Application Interface

The first page the user sees when logged in to the mobile application is a list of active Device IDs. An active device has been initiated (i.e. turned on) in the last 24 hours. To refresh the list of active devices, swipe down.

Settings

Close Settings	
🗘 General	
Permissions >	
Tutorial >	
E Legal	
Legal Notices	
Implement (01) 07540162030048 (10) Implement MedEnvoy Prinses Margrietplantsoen 33, Suite 123 2595 AM The Hague, The Netherlands io.my01.ios-client@ /Production/Release	

Instructions

Connecting to the MY01 device

Open the MY01 Mobile Application and pull down to scan for nearby advertising MY01 devices. Each entry in the list contains the unique MY01 Device identifier.

Verify that the unique identifier displayed on the application matches the one displayed on the MY01 device's LCD top-left corner.

Tap on a device entry to connect to the MY01 device. If connecting for the first time, a Bluetooth pairing request will appear. Confirm that the pairing PIN displayed on the application matches the pairing PIN displayed on the MY01 device's LCD. Click on the "Pair" button to confirm and connect.

Once connected to a MY01 device, a Bluetooth mark will appear next to the MY01 device unique identifier. Pressure data transmission will begin shortly after connection.

(10:40		0	add ute 📼	١
	=	MY01 In	ic. ~	P	
1		nt data is not f patient monito	or diagnos oring.		
Ш	• A1A1 *	Bluetooth Mark	mmHg	7.5	
	Streaming.				
I					I
	_		2		

Notes :

- The application will try to automatically reconnect to nearby advertising devices. In certain cases the application won't be able to connect automatically. For more information on what to do, refer to the troubleshooting section.
- Only one phone can be connected to a device at a time.
- If the application doesn't connect to the device, refer to the troubleshooting section for help.

Active Device information

Data Display

Warning : do NOT use the MY01 device or the MY01 Mobile Application for remote monitoring or for diagnosis.

Important: in order for the recorded data to be successfully transferred to the Cloud for later retrieval, the user MUST ensure that The MY01 device connected the the MY01 Mobile Application.

Once the MY01 device has been successfully connected to the MY01 Mobile Application the graphing function can be accessed by pressing on the "Device of choice" section of the interface. A Bluetooth symbol will be displayed if the device is connected and streaming data.

10:40		⇒ ∘	add ure 📼
=	M	/01 Inc. 🕚	~ R
Displayed	patient data is n patient mo	ot for diagnosis initoring.	or active
READY TO	O CONNECT		
• B2B2		mmHg	
MAY 5 20	20		
• E9X7	*	mmHg 2	2.1
Streami	ng.	دالس	
		\cup	
	.	2	

To change the display from portrait to landscape mode, simply rotate your mobile device by 90 degrees and the screen will automatically rotate.

Switching Organizations

Users that are registered with multiple organizations are able to toggle between different organizations by clicking on the drop down list at the top of the home screen.

Note: the drop down arrow will only appear for users who are registered in more than one organization.

Note: a device can only be registered within one organization. Should you try to connect to a device that has been previously registered to another organization the warning below will appear. Ensure that you have selected the correct organization and connect to the device again.

Disconnecting from a MY01 device

To disconnect from a MY01 device, slide the device entry to the left and click on the "Disconnect" button. The Bluetooth logo will disappear from the MY01 Device.

List of inactive devices

A list of the inactive devices that have been used in the past is shown and pressure measurements can be viewed.

Logging out of the MY01 Mobile Application

Once users have finished using the MY01 device, they can log out of the mobile application by tapping opening the "Menu" and clicking on "Log Out".

Help & Troubleshooting

Help Center/ FAQs

We have compiled a list of tutorials and documentation to help users and respond to any questions that should arise. The "Help Center" can be accessed from the Menu screen, located on the upper left of the Home screen.

Chat

Should you not find the required information on our help center we invite you to go to our chat system. A MY01 Representative can be contacted directly in the MY01 Mobile Application by pressing the "Chat" button on the upper right of the home screen.

10:40	— •			10:40	— •	II LTE 🗩
≡ M	Y01 Inc. 🗸	P.	<		MY01 Inc.	×
Displayed patient data patien	is not for diagnos t monitoring.	sis or active		Ask us a feedback. V	anything, or share you Ve're happy to help ge I on what we're workin	ur et you
READY TO CONNEC	т		h 1			
• B2B2	mmHg					
MAY 5 2020				Will rep	oly as soon as they can	
• A1A1 *	mmHg	0.2				
Streaming.						
			Sta	art a conversati	on	
±	2	a	Aa	GIF 💭		Send
		•				

Issue	Resolution
Application will not start	Delete and reinstall the app.
	 Check for application updates in the applicable digital distribution platform (Application store & Play Store).
Application crashes	Close the application and restart the app.
	Check for application updates in the applicable digital distribution platform (Application store & Play Store).
Application cannot connect to a device Device does not appear in the device list, or does not appear in the "Ready to Connect" section	 Ensure the device is turned on and in the same room, within 4 meters of the mobile device. Ensure the device's LCD is not displaying any warning. Ensure your mobile device has Bluetooth turned on. Ensure your mobile device has network connectivity Active / ON (WiFI or cellular connection). Ensure the device is not already connected to another app. Ensure that the Bluetooth logo is not displayed on the device's LCD Refresh the device list by "pulling down". If the steps above are not successful. Close and restart the app. Check for application updates in the applicable application store.
Application cannot connect to a device. When a device entry is pressed, the application displays the device's pressure chart instead of connecting via Bluetooth.	 Ensure the device entry appears in the "Ready to Connect" section before pressing it. If the device entry is not in the "Ready to Connect" section, refer to "Device does not appear in the device list, or does not appear in the "Ready to Connect" section" above.
Application cannot connect to a device An error is displayed upon connection	 Ensure the device is turned on and in the same room, within 4 meters of the mobile device. Ensure the device's LCD is not displaying any warning. Ensure your mobile device has network connectivity turned ON (WiFI or cellular connection). Open your mobile device's Bluetooth settings and "forget" the appropriate device. Connect by pressing the device entry under the "Ready To Connect" section of the active devices list. If the steps above are not successful, close and restart the app. Check for application updates in the applicable digital distribution platform (Application store & Play Store).
A device with incompatible firmware has been scanned	 Ensure you are using the latest version of MY01 Mobile Application. Try to reconnect to the device
Cannot find a device turned on more than 24 hours ago.	 Ensure your mobile device has network connectivity (WiFI or cellular connection). Ensure the inactive devices list is currently displayed. Press the left

	 tab on the main application screen. Refresh the device list by "pulling down". If not successful, it may indicate that MY01 service is undergoing maintenance or is unavailable. Please try again later.
Application displays "Bluetooth is off."	 Turn on bluetooth on your mobile device. If not successful, close the application and restart it.
Application displays "Bluetooth is not supported on this mobile device."	 Ensure you are using a Bluetooth 4.2 compatible mobile device. For the latest list of compatible mobile device open the in-application help center or visit <u>https://help.my01.io</u>
Application opens but cannot login	 Ensure your mobile device has network connectivity (WiFI or cellular connection). Try to re-enter your credentials. You can use the "eye" icon to reveal your password before submitting. Press the "Forgot Password?" button on the login screen or open https://gateway.my01.io/account/reset/request using a web browser.
Application displays "Mobile device not secure" when started.	 Indicates that your mobile device is not protected by a passcode. You can configure a passcode in your mobile device's security settings.
Application displays "Server unreachable"	 Ensure your mobile device has network connectivity (WiFI or cellular connection). If not successful, it may indicate that MY01 service is undergoing maintenance or is unavailable. Please try again later.

Technical Issues

Paper IFU manuals are available upon request free of charge, within 7 days of the request. If you want to request a paper copy of the IFU, please email help@MY01.io or call +1 (855) 799-6901.

For bug reporting, technical difficulties, any related quality issues, or any serious incident that has occurred in relation to the device contact: please contact help@MY01.io or call +1 (855) 292-6901.

For additional resources, visit <u>help.my01.io</u> or call +1 (855) 799-6901

A MY01 representative will deal with any quality issues related to hardware or software functionality in a timely manner.

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