
# User Manual / Guide d'utilisation MY01 Mobile Application/ Application mobile MY01

MY01 Inc.  
400 Boul de Maisonneuve Ouest,  
Suite 700, Montréal, QC, Canada, H3A 1L4  
+1 (855) 799-6901

MedEnvoy  
Prinses Margrietplantsoen 33  
Suite 123, 2595 AM The Hague  
The Netherlands

Rx ONLY

CEpartner4U BV  
Esdoornlaan 13  
3951 DB Maarn  
The Netherlands

## Table of Contents

<b>Introduction</b> .....	<b>5</b>
<b>Getting Started</b> .....	<b>6</b>
Recommended smart device settings.....	6
Bluetooth.....	6
Date & Time .....	7
Signing up .....	8
Signing in .....	8
Application Interface.....	9
<b>Settings</b> .....	<b>10</b>
<b>Instructions</b> .....	<b>11</b>
Connecting to the MY01 device .....	11
Active Device information.....	13
Data Display.....	14
Switching Organizations.....	15
Disconnecting from a MY01 device .....	17
List of inactive devices .....	18
Logging out of the MY01 Mobile Application.....	19
<b>Help &amp; Troubleshooting</b> .....	<b>20</b>
Help Center/ FAQs.....	20
Chat .....	21
Technical Issues.....	23
<b>Présentation</b> .....	<b>26</b>

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Symbol Explanation	Symbol
Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	
Manufacturer	
Authorized representative in the European Community	
Medical Device	
Unique Device Identifier (UDI)	
Importer	

# Introduction

## **Intended use**

The MY01 Mobile Application is an optional application intended for storing and displaying identical pressure values from the MY01 Continuous Compartmental Pressure Monitor device. The data is for informational purposes only and is not intended to be used for diagnosis of any nature or active patient monitoring.

## **Patient target group**

The device is intended for use in patients at risk or with suspected compartment syndrome due to a wide variety of clinical presentations including trauma.

## **Intended Users**

Trained trauma surgeons/Physicians.

## **Application download**

# Getting Started

## Recommended smart device settings

Ensure your smart device settings correspond with the following parameters:

**Compatibility:** For a list of smart devices and operating systems that work with the MY01 Mobile Application check [help.my01.io](http://help.my01.io)

- Before downloading and installing the app, make sure your mobile device meets or exceeds the minimum performance specifications. Additional information on the most up to date device requirements can be found at [help.my01.io](http://help.my01.io)
- The MY01 device uses Bluetooth Low Energy (BLE) technology specification 4.2 Ensure that the mobile phone uses BLE-compatible wireless technology.
- Recommended: Access to your mobile device should be password protected.

## Bluetooth

Enabling Bluetooth on your iOS device

1. Tap on **Settings** on your iPhone.
2. Select **General**.
3. Tap on **Bluetooth**.
4. Toggle Bluetooth to **On**.

Enabling Bluetooth on your Android device:

1. Tap on **Settings** on your Android device.
2. Look for **Bluetooth** or the Bluetooth symbol in your settings and tap it.
3. There should be an option to enable. Please tap or swipe on it so that it is in the **on** position.

## Date & Time

Update Date & Time on Your iOS Device

1. From your home screen, navigate to Settings.
2. Tap General.
3. Tap Date & Time.
4. Make sure that the option Set Automatically is turned on.
5. If this option is turned off, check that the correct Date, Time and Time Zone are selected.

Update Date & Time on Your Android Device

1. Tap Settings to open the Settings menu.
2. Tap Date & Time.
3. Tap Automatic.
4. If this option is turned off, check that the correct Date, Time and Time Zone are selected.

## Signing up

To connect the MY01 device through Bluetooth to the MY01 Mobile Application for the first time, users must follow these steps:

1. Users are invited by their institution to join the MY01 Mobile Application. Choose a secure password and confirm your account.
2. Install the MY01 Mobile Application from the Google Play Store or Apple App Store.

## Signing in

1. Open the MY01 Mobile Application.
2. Log in using your credentials.
3. If you've forgotten your password, click on reset password and follow the on-screen instructions.

# Application Interface

The first page the user sees when logged in to the mobile application is a list of active Device IDs. An active device has been initiated (i.e. turned on) in the last 24 hours. To refresh the list of active devices, swipe down.

**Menu**

**Help**  
Connects you to a MY01 representative



**Ready to Connect**  
Lists devices which are turned on and ready to pair with the mobile device

**In Use**  
Lists devices that are currently transmitting data to the mobile application.

**Active / Inactive Device Tabs**  
The user can navigate between lists of active and inactive devices by selecting their respective tabs.

10:40  
MY01 Inc. v  
Displayed patient data is not for diagnosis or active patient monitoring.  
**READY TO CONNECT**  
• A1A1 mmHg --  
MAY 5 2020  
• A1A1 \* mmHg 0.2  
Streaming.  
[User Icon] [Search Icon]

# Settings



# Instructions


## Connecting to the MY01 device

Open the MY01 Mobile Application and pull down to scan for nearby advertising MY01 devices. Each entry in the list contains the unique MY01 Device identifier.


Verify that the unique identifier displayed on the application matches the one displayed on the MY01 device's LCD top-left corner.

Tap on a device entry to connect to the MY01 device. If connecting for the first time, a Bluetooth pairing request will appear. Confirm that the pairing PIN displayed on the application matches the pairing PIN displayed on the MY01 device's LCD. Click on the "Pair" button to confirm and connect.


Once connected to a MY01 device, a Bluetooth mark will appear next to the MY01 device unique identifier. Pressure data transmission will begin shortly after connection.


## Notes:

- The application will try to automatically reconnect to nearby advertising devices. In certain cases, the application won't be able to connect automatically. For more information on what to do, refer to the troubleshooting section.
- Only one phone can be connected to a device at a time.
- If the application doesn't connect to the device, refer to the troubleshooting section for help.

## Active Device information


## Data Display


**Warning:** Do not use the MY01 Mobile Application for remote monitoring or for diagnosis.

**Important:** In order for the recorded data to be successfully transferred to the Cloud for later retrieval, the user must ensure that the MY01 device is connected to the MY01 Mobile Application.

Once the MY01 device has been successfully connected to the MY01 Mobile Application the graphing function can be accessed by pressing on the “Device of choice” section of the interface. A Bluetooth symbol will be displayed if the device is connected and streaming data.


To change the display from portrait to landscape mode, simply rotate your mobile device by 90 degrees and the screen will automatically rotate.


## Switching Organizations

Users that are registered with multiple organizations are able to toggle between different organizations by clicking on the drop-down list at the top of the home screen.

**Note:** the drop-down arrow will only appear for users who are registered in more than one organization.


Note: a device can only be registered within one organization. Should you try to connect to a device that has been previously registered to another organization the warning below will appear. Ensure that you have selected the correct organization and connect to the device again.


## Disconnecting from a MY01 device

To disconnect from a MY01 device, slide the device entry to the left and click on the “Disconnect” button. The Bluetooth logo will disappear from the MY01 Device.


# List of inactive devices

A list of the inactive devices that have been used in the past is shown and pressure measurements can be viewed.

## Logging out of the MY01 Mobile Application



Once users have finished using the MY01 device, they can log out of the mobile application by tapping opening the “Menu” and clicking on “Log Out”.

# Help & Troubleshooting


## Help Center/ FAQs

We have compiled a list of tutorials and documentation to help users and respond to any questions that should arise. The “Help Center” can be accessed from the Menu screen, located on the upper left of the home screen.

# Chat

Should you not find the required information on our help center we invite you to go to our chat system. A MY01 Representative can be contacted directly in the MY01 Mobile Application by pressing the “Chat” button on the upper right of the home screen.

Issue	Resolution
Application will not start	<ul style="list-style-type: none"> <li>● Delete and reinstall the app.</li> <li>● Check for application updates in the applicable digital distribution platform (Application store &amp; Play Store).</li> </ul>
Application crashes	<ul style="list-style-type: none"> <li>● Close the application and restart the app.</li> <li>● Check for application updates in the applicable digital distribution platform (Application store &amp; Play Store).</li> </ul>
<p>Application cannot connect to a device</p> <p>Device does not appear in the device list, or does not appear in the “Ready to Connect” section</p>	<ul style="list-style-type: none"> <li>● Ensure the device is turned on and in the same room, within 4 meters of the mobile device.</li> <li>● Ensure the device’s LCD is not displaying any warning.</li> <li>● Ensure your mobile device has Bluetooth turned on.</li> <li>● Ensure your mobile device has network connectivity Active / ON (WiFi or cellular connection).</li> <li>● Ensure the device is not already connected to another app. Ensure that the Bluetooth logo is not displayed on the device’s LCD</li> <li>● Refresh the device list by “pulling down”.</li> <li>● If the steps above are not successful. Close and restart the app. Check for application updates in the applicable application store.</li> </ul>
Application cannot connect to a device. When a device entry is pressed, the application displays the device’s pressure chart instead of connecting via Bluetooth.	<ul style="list-style-type: none"> <li>● Ensure the device entry appears in the “Ready to Connect” section before pressing it.</li> <li>● If the device entry is not in the “Ready to Connect” section, refer to “Device does not appear in the device list, or does not appear in the “Ready to Connect” section” above.</li> </ul>
<p>Application cannot connect to a device</p> <p>An error is displayed upon connection</p>	<ul style="list-style-type: none"> <li>● Ensure the device is turned on and in the same room, within 4 meters of the mobile device.</li> <li>● Ensure the device’s LCD is not displaying any warning.</li> <li>● Ensure your mobile device has network connectivity turned ON (WiFi or cellular connection).</li> <li>● Open your mobile device’s Bluetooth settings and” forget” the appropriate device.</li> <li>● Connect by pressing the device entry under the “Ready to Connect” section of the active devices list.</li> <li>● If the steps above are not successful, close and restart the app.</li> <li>● Check for application updates in the applicable digital distribution platform (Application store &amp; Play Store).</li> </ul>
A device with incompatible firmware has been scanned	<ul style="list-style-type: none"> <li>● Ensure you are using the latest version of MY01 Mobile Application.</li> <li>● Try to reconnect to the device</li> </ul>

Cannot find a device turned on more than 24 hours ago.	<ul style="list-style-type: none"> <li>• Ensure your mobile device has network connectivity (Wi-Fi or cellular connection).</li> <li>• Ensure the inactive devices list is currently displayed. Press the left tab on the main application screen.</li> <li>• Refresh the device list by “pulling down”.</li> <li>• If not successful, it may indicate that MY01 service is undergoing maintenance or is unavailable. Please try again later.</li> </ul>
Application displays “Bluetooth is off.”	<ul style="list-style-type: none"> <li>• Turn on Bluetooth on your mobile device. If not successful, close the application and restart it.</li> </ul>
Application displays “Bluetooth is not supported on this mobile device.”	<ul style="list-style-type: none"> <li>• Ensure you are using a Bluetooth 4.2 compatible mobile device. For the latest list of compatible mobile device open the in-application help center or visit <a href="https://help.my01.io">https://help.my01.io</a></li> </ul>
Application opens but cannot login	<ul style="list-style-type: none"> <li>• Ensure your mobile device has network connectivity (Wi-Fi or cellular connection).</li> <li>• Try to re-enter your credentials. You can use the “eye” icon to reveal your password before submitting.</li> <li>• Press the “Forgot Password?” button on the login screen or open <a href="https://gateway.my01.io/account/reset/request">https://gateway.my01.io/account/reset/request</a> using a web browser.</li> </ul>
Application displays “Mobile device not secure” when started.	<ul style="list-style-type: none"> <li>• Indicates that your mobile device is not protected by a passcode. You can configure a passcode in your mobile device’s security settings.</li> </ul>
Application displays “Server unreachable”	<ul style="list-style-type: none"> <li>• Ensure your mobile device has network connectivity (Wi-Fi or cellular connection).</li> <li>• If not successful, it may indicate that MY01 service is undergoing maintenance or is unavailable. Please try again later.</li> </ul>

## Technical Issues

Paper IFU manuals are available upon request free of charge, within 7 days of the request. If you want to request a paper copy of the IFU, please email [info@MY01.io](mailto:info@MY01.io) or call +1 (855) 799-6901.

For bug reporting, technical difficulties, any related quality issues, or any serious incident that has occurred in relation to the device contact: please contact [info@MY01.io](mailto:info@MY01.io) or call +1 (855) 799-6901.

For additional resources, visit [help.my01.io](https://help.my01.io) or call +1 (855) 799-6901. A MY01 representative will deal with any quality issues related to hardware or software functionality in a timely manner.

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