# User Manual

---

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Continuous Compartmental Pressure Monitor

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This device is protected by patents and has additional patents pending.

**Device Name:** MY01 Continuous Compartmental Pressure Monitor

**Model number:** MY01-0001

**R<sub>x</sub> ONLY**

MedEnvoy  
Prinses Margriplantsoen 33  
Suite 123, 2595 AM The Hague  
The Netherlands  
[info@medenvoyglobal.com](mailto:info@medenvoyglobal.com)  
[medenvoyglobal.com](http://medenvoyglobal.com)

MY01 Inc.  
400 Boul de Maisonneuve Ouest,  
Suite 700, Montréal, QC H3A 1L4  
+1 (855) 799-6901  
[info@my01.io](mailto:info@my01.io) / [my01.io](http://my01.io)

MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland  
[info@medenvoyglobal.com](mailto:info@medenvoyglobal.com)  
[medenvoyglobal.com](http://medenvoyglobal.com)

CEpartner4U BV  
Esdoormlaan 13  
3951 DB Maarn  
The Netherlands  
[cepartner4u.com](http://cepartner4u.com)  
[office@cepartner4u.com](mailto:office@cepartner4u.com)

# TABLE OF CONTENTS

|                                                    |           |
|----------------------------------------------------|-----------|
| <b>Glossary of Labeling Symbols .....</b>          | <b>4</b>  |
| Device Labeling Symbols .....                      | 4         |
| Device Interface Symbols .....                     | 5         |
| <b>Introduction.....</b>                           | <b>6</b>  |
| Warnings and Notes.....                            | 6         |
| Intended Use.....                                  | 6         |
| <b>User &amp; Patient Safety Precautions .....</b> | <b>8</b>  |
| <b>Components.....</b>                             | <b>9</b>  |
| <b>Device Interface.....</b>                       | <b>10</b> |
| <b>Wireless Connection .....</b>                   | <b>12</b> |
| <b>Operational Instructions.....</b>               | <b>14</b> |
| <b>Troubleshooting .....</b>                       | <b>21</b> |
| Device Interface Error Codes .....                 | 21        |
| Instrument Defects or Product Deficiencies .....   | 22        |
| <b>Specifications .....</b>                        | <b>23</b> |
| Device .....                                       | 23        |
| Environmental Conditions .....                     | 24        |
| Electromagnetic Compatibility .....                | 25        |
| <b>Reorder Information.....</b>                    | <b>28</b> |

# Glossary of Labeling Symbols

## DEVICE LABELING SYMBOLS

| Symbol                                                                              | Title and Description                                                                                                                                                                                                                                                                                                | Reference           |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
|    | Caution                                                                                                                                                                                                                                                                                                              | ISO 15223-1 (5.4.4) |
|    | Follow instructions for use                                                                                                                                                                                                                                                                                          | ISO 7010-M002       |
| <b>Rx ONLY</b>                                                                      | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician                                                                                                                                                                                                                           | 21 CFR 801.109      |
|     | Sterilized using ethylene oxide                                                                                                                                                                                                                                                                                      | ISO 15223-1 (5.2.3) |
|    | Type BF applied part                                                                                                                                                                                                                                                                                                 | IEC 60417 (5333)    |
|    | Do not resterilize                                                                                                                                                                                                                                                                                                   | ISO 15223-1 (5.2.6) |
|   | Do not use if package is damaged                                                                                                                                                                                                                                                                                     | ISO 15223-1 (5.2.8) |
|  | Do not re-use. Single use only                                                                                                                                                                                                                                                                                       | ISO 15223-1 (5.4.2) |
|  | Manufacturer                                                                                                                                                                                                                                                                                                         | ISO 15223-1 (5.1.1) |
|  | Maximum and minimum temperature limits                                                                                                                                                                                                                                                                               | ISO 15223-1 (5.3.7) |
|  | Maximum and minimum relative humidity limits                                                                                                                                                                                                                                                                         | ISO 15223-1 (5.3.8) |
|  | Maximum and minimum pressure limits                                                                                                                                                                                                                                                                                  | ISO 15223-1 (5.3.9) |
|  | Keep away from sunlight                                                                                                                                                                                                                                                                                              | ISO 15223-1 (5.3.2) |
|  | Use-by date                                                                                                                                                                                                                                                                                                          | ISO 15223-1 (5.1.4) |
| <b>LOT</b>                                                                          | Batch code                                                                                                                                                                                                                                                                                                           | ISO 15223-1 (5.1.5) |
| <b>REF</b>                                                                          | Catalogue number                                                                                                                                                                                                                                                                                                     | ISO 15223-1 (5.1.6) |
| IPn1n2                                                                              | Ingress Protection (IP Rating):<br>The first digit indicates the level of protection that the enclosure provides against the ingress of solid foreign objects (Scale: 0-6).<br>The second digit indicates the level of protection of the equipment inside the enclosure against the ingress of liquids (Scale: 0-8). | IEC 60529           |

|  |                                                           |                                       |
|--|-----------------------------------------------------------|---------------------------------------|
|  | Consult instructions for use                              | ISO 15223-1 (5.4.3)                   |
|  | Serial number                                             | ISO 15223-1 (5.1.7)                   |
|  | Unique device identifier                                  | ISO 15223-1 (5.7.10)                  |
|  | Medical device                                            | ISO 15223-1 (5.7.7)                   |
|  | Model Number                                              | ISO 15223-1 (5.1.10)                  |
|  | Signifies conformity with European legislation            | Regulation (EU) 2017/745 (Article 20) |
|  | Keep away from magnetic resonance imaging (MRI) equipment | ASTM F2503-23                         |
|  | Importer                                                  | ISO 15223-1 (5.1.8)                   |
|  | Indicates the authorized representative in Switzerland    | ISO 15223-1 (5.1.2)                   |
|  | Authorized representative in the European Community       | ISO 15223-1 (5.1.2)                   |

## DEVICE INTERFACE SYMBOLS

| Symbol | Title and Description                                     | Reference                                                                                                       |
|--------|-----------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
|        | Bluetooth                                                 | Trademark of Bluetooth Special Interest Group (Bluetooth SIG, Inc.), Bluetooth ® wireless or enabled technology |
|        | Full battery                                              | Proprietary                                                                                                     |
|        | Warning: Low battery level                                | Proprietary                                                                                                     |
|        | Trend arrows: indicates rising or falling pressure values | Proprietary                                                                                                     |

# Introduction

## WARNINGS AND NOTES

This User Manual is intended for medical professionals who have completed appropriate MY01 training within their respective healthcare facilities. It provides information to ensure safe and effective use of the product. The following terminology is used:

**WARNING: The personal safety of the patient may be involved. Disregarding this information could result in injury to the patients, users, and/or others.**

**CAUTION:** This alerts the user to a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the patient, user, or others, or cause damage to the device or other property. It may also indicate the need for special care for safe and effective use of the device.

**NOTE:** This provides additional important information the user should be aware of.

## INTENDED USE

### CANADA AND UNITED-STATES

The MY01 Continuous Compartmental Pressure Monitor is used for real-time and continuous measurement of the muscle compartment pressure. The measured muscle compartment pressure can be used as an aid in diagnosis of Compartment Syndrome (Acute and Chronic). The MY01 Mobile Application is an application intended for storing and displaying identical pressure values from the MY01 Continuous Compartmental Pressure Monitor device and calculating critical muscle perfusion pressure utilizing diastolic pressure manual entry by the physician. Diagnosis should always be made in conjunction with clinical assessments.

**NOTE:**

- Refer to MY01 Mobile Application User Manual (REF: MYO-00566-M)
- The MY01 App and MY01 device must be used in tandem.

## IN CONFORMITY WITH EUROPEAN LEGISLATION

The MY01 Continuous Compartmental Pressure Monitor is intended for real-time and continuous measurement of compartmental pressures. The measured compartmental pressures can be used as an aid in the diagnosis of compartment syndrome. The trend arrows displayed are meant for qualitative purposes only and are not intended to have any clinical significance. The MY01 Mobile Application is an optional application intended for storing and displaying identical pressure values from the MY01 Continuous Compartmental Pressure Monitor device. The data is for informational purposes only and is not intended to be used for diagnosis of any nature or active patient monitoring.

NOTE: The MY01 Mobile Application that is currently available in conformity with European Legislation does not include functionality for the display of muscle perfusion pressure; see MY01 Mobile Application User Manual (REF: MYO-00319-M)

**R<sub>x</sub> ONLY**

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

# User & Patient Safety Precautions

## WARNING

- The device should only be used in a healthcare facility environment.
- Users (Medical Professionals) must complete the required training **BEFORE** attempting to use the device to ensure safe and effective use. Training package can be accessed here <https://qrco.de/my01-training> or by scanning the QR code.
- Do not perform a diagnosis solely based on pressure measurements of the device. Always use the device with the current standard of care.
- Use aseptic practices during usage - follow healthcare facility guidelines.
- Do not ship, store or use the device outside the specified environmental conditions (See Table 4).
- Single-use device, do not reuse. Patient safety may be compromised.
- Do not use the device past the expiration date.
- Do not resterilize the device. It is initially ETO sterilized and cannot be reprocessed.
- Do not clean the device.
- Do not perform any maintenance on the device.
- Do not use the same device in multiple patients.
- Do not use the device in proximity to MRI equipment and high frequency surgical equipment.
- Do not use the device if it is believed to be faulty. The user should exercise clinical judgment when performing measurements.

## CONTRAINDICATIONS

No known contraindications

# Components

## Packaging and Device

Figure 1: Packaged MY01 Device

The MY01 Continuous Compartmental Pressure Monitor comprises the following 2 major components: the Introducer and the Pressure Monitor.


Figure 2: Introducer


Figure 3: Pressure Monitor (Applied Part)

# Device Interface

Figure 4: Device Controls and Display Icons

## Significance of trend arrows

The decrease in pressure is greater than or equal to 0.5 mmHg/hour.

The decrease in pressure is greater than or equal to 2 mmHg/hour.

The increase in pressure is greater than or equal to 0.5 mmHg/hour.

The increase in pressure is greater than or equal to 2 mmHg/hour.

 **WARNING**

- Trend arrows are meant for qualitative purposes only and are not clinically significant. Do not use trend arrows for diagnostic purposes. Always use pressure measurements and clinical judgment along with the current standard of care.

Figure 5: Rotating Display Orientation

# Wireless Connection

## BLUETOOTH SPECIFICATIONS

The MY01 device transmits identical pressure readings collected over time to the MY01 Mobile Application using an authenticated BLE link. The MY01 device uses Bluetooth Low Energy (BLE) technology specification 4.2. Ensure that the mobile phone uses BLE-compatible wireless technology.

Table 1: BLE Specifications

|                                     |                                                                                      |
|-------------------------------------|--------------------------------------------------------------------------------------|
| <b>Wireless technology version</b>  | BLE specification 4.2                                                                |
| <b>Frequency</b>                    | 2.4GHz                                                                               |
| <b>RF frequencies bandwidth</b>     | 2.402GHz - 2.4835GHz                                                                 |
| <b>RF radiated output power</b>     | 0dBm, 1mWatts                                                                        |
| <b>Modulation type</b>              | Gaussian Frequency Shift Keying (GFSK)                                               |
| <b>Bluetooth Interface Security</b> | Authenticated LE Secure Connections pairing uses Numeric Comparison with P-256 ECDH. |
|                                     | Application layer Device-Application authentication using AES encryption             |

## MY01 MOBILE APPLICATION

The MY01 Mobile Application (MY01 App) provides a graphical user interface for storing and displaying identical pressure values from the MY01 device. The version of the MY01 App that is authorized for use in the US and Canada, in addition to the muscle compartment pressure, calculates and displays the muscle perfusion pressure over time based on manual entries of the patient's diastolic blood pressure. Diagnosis should always be made in conjunction with clinical assessments.

For directions on how to use the MY01 App, please refer to the MY01 Mobile Application User Manual that is available in your region:

- MY01 App with compartmental and perfusion pressure display (P/N: MYO-00566-M);
- MY01 App with compartmental pressure display only (P/N: MYO-00319-M).

**NOTE:**

- The MY01 App does not control the MY01 device in any way.

## MY01 MOBILE APPLICATION DOWNLOAD

# Operational Instructions

**Step 1 :** Peel the Tyvek lid where indicated. Lift the top cover and remove the device from the packaging. While holding the Introducer (see Figure 6), remove the Needle Cap before turning on the device.

## WARNING

- Do not use the device if the Sensor is found unhooked from the Needle after opening the package. Never attempt to re-assemble the device.
- Do not use the device if the integrity of the packaging is compromised.

**Step 2 :** Activate the Pressure Monitor by pressing and holding the Push-Button until the MY01 logo appears. The Display will subsequently show the MY01 logo and a pressure value of 0.

## NOTE

- Make sure to not touch the sensor tip while turning ON the Device.

Figure 6: Removing Needle Cap and Turning ON the Pressure Monitor

## WARNING

- Always turn ON the device BEFORE introducing into patient.
- Do not use the device if pressure measurements are not within -2.5 to 2.5 mmHg prior to insertion.


Figure 7: Recommended Device Grips

**Step 3** : Without applying excessive force on the Display, insert the Needle gently into the muscle compartment in a controlled linear motion. Markings on the Needle should be used to estimate the depth of the insertion.


Figure 8: Introducing the Device

## CAUTION

- Prior to insertion, identify the target anatomical compartment based on clinical assessment. The selected compartment must provide sufficient depth and volume to accommodate the MY01 needle's length and diameter. The sensor must be positioned within the intended compartment to ensure reliable intracompartmental pressure measurements. Improper placement may result in inaccurate readings.

#### NOTE

- The Display can be rotated to the desired orientation to ensure visibility during the insertion (see Figure 5).
- Single marks are spaced along the needle at 1 cm intervals with the first located 2 cm from the Needle tip. Double marks are placed at 5 cm intervals as measured from the Needle tip.
- Consider applying local anesthetic before insertion. Ensure there are no allergies to anesthetic use.
- If necessary, a small incision can be made to facilitate the insertion process.
- The Introducer can be retracted and re-inserted (up to 5 total insertions) to perform additional single-point measurements without ejecting the Pressure Monitor.

#### WARNING

- Always keep the Atmospheric Pressure Vents on the bottom-right of the screen unobstructed (see Figure 4).
- Do not use the device if the Sensor unhooks unexpectedly after removing the Needle Cap or at any point during the insertion (up to 5 single-point insertions). Never attempt to re-assemble the device.
- Do not rotate the Introducer during insertion to prevent premature Sensor unhooking.

**Step 4 :** When the Sensor is in the desired position and readings have stabilized, eject the Pressure Monitor from the Introducer by pressing gently through the back opening of the Introducer. An adhesive strip will be exposed on the back of the Device Body when ejected.

A

B

Figure 9: Ejecting the Pressure Monitor

**Step 5 :** While holding the Introducer in one hand, use the other hand to adhere the Device Body to the patient's skin using the exposed adhesive strip on the back of the Device Body. Position the Device Body face-up on the patient's skin near the insertion point, ensuring sufficient slack in the Lead-Wire. The Lead-Wire should extend straight out from the insertion site indicating the insertion angle of the Sensor.


Figure 10: Placing the Device Body on Patient Skin

#### NOTE

- Care should be taken when moving the Device Body to not pull the Lead-Wire which can displace the Sensor within the muscle compartment.
- The Device Body should be installed in a location which won't interfere with adhesion during the monitoring period.
- It is recommended to prepare the skin appropriately to improve adhesion, especially on hairy or oily skin.

#### WARNING

- If the adhesive strip is not exposed after ejecting the Pressure Monitor from the Introducer (see Figure 9), manually remove the protective liner. Do not re-apply the adhesive. Additional medical tape should be used if the location of the Device Body is changed after the first application. Ensure that the Atmospheric Pressure Vents (see Figure 4) are not covered by the medical tape.

**Step 6 :** To disengage the Sensor, use a finger to firmly press the lead-wire against the skin on the insertion site, then rotate the Introducer 180 degrees. While keeping the finger pressed on the lead-wire, gently remove the Introducer from the patient. After removing the Introducer, dispose of it in a biohazard-sharps receptacle, following facility guidelines.

#### NOTE

- It is recommended to apply a dressing to the insertion site to fix the Lead-Wire in place.

#### WARNING

- Sharps biohazard - dispose as per facility guidelines and/or local regulations.


Figure 11: Removing the Introducer

**Step 7** : Monitor pressure readings for a period of up to 18 hours. The Pressure Monitor should be routinely checked to ensure that it is secured to the patient throughout the monitoring period and that the Lead-Wire does not pull on the Sensor and displace it.

 **WARNING**

- The Pressure Monitor may stop functioning correctly after 18 hours when the battery is depleted. Do not leave the Sensor inside the patient for a period longer than 18 hours.


**Step 8** : When monitoring is complete, remove the dressing and gently pull on the Lead-Wire by hand to remove the Sensor from the patient. Dispose the Pressure Monitor in a biohazard container, as per facility guidelines.

**NOTE**

- The Lead-Wire should be pulled out at the same angle used for insertion of the Sensor to minimize removal forces.

 **WARNING**

- The device is for single-use only (a single monitoring session). Do not attempt to reassemble the device, to reuse or replace the batteries after use.
- Take necessary precautions while using the device intraoperatively (e.g. fasciotomy, amputation, debridement, skin graft).

A


B

C

D

E


F

G

H

I

Figure 12: Insertion Procedure Summary

# Troubleshooting

## DEVICE INTERFACE ERROR CODES

Table 2: Error Code List

| ERROR CODE          | ERROR NAME                       | DESCRIPTION                                                                    | ACTION                                                                                                               |
|---------------------|----------------------------------|--------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| NON-CRITICAL ERRORS |                                  |                                                                                |                                                                                                                      |
| 3                   | Lifetime Exceeded                | Device has been running for longer than intended operating time                | Device operation is not recommended beyond this time, if device is currently being used cease operations immediately |
| 4                   | Outside Operating Pressure Range | The absolute pressure reading from sensors is outside operating pressure range | Disregard the pressure values while error sign is visible                                                            |
| CRITICAL ERRORS     |                                  |                                                                                |                                                                                                                      |
| ERR 10              | Sensor Integrity Compromised     | System has detected a critical malfunction in the sensing element              | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 11              | Software Error                   | Software malfunction                                                           | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 101             | Maximum Zeroing Entries          | May occur if sensor tip is touched during activation.                          | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 201 - 204       | Sensor Tip Errors                | Sensor experiencing internal electrical or signal issues.                      | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 205             | Sensor Communication Lost        | Device is unable to communicate with sensor tip hardware.                      | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 301             | Display Unresponsive             | LCD interface failed to initialize or respond                                  | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 401             | System Overheating               | Internal chip temperature is above the safe threshold for operation            | DO NOT use the device. Call customer support for instructions on how to return the device.                           |
| ERR 402             | System Too Cold                  | Device operation environment is too cold                                       | DO NOT use the device. Call customer support for instructions on how to return the device.                           |

## CRITICAL ERRORS

|         |                            |                                                           |                                                                                            |
|---------|----------------------------|-----------------------------------------------------------|--------------------------------------------------------------------------------------------|
| ERR 403 | Memory Corruption Detected | Device memory may have been interrupted during operation. | DO NOT use the device. Call customer support for instructions on how to return the device. |
| ERR 404 | System Resets              | Device has reset multiple times unexpectedly.             | DO NOT use the device. Call customer support for instructions on how to return the device. |
| ERR 405 | Invalid System State       | Device attempted an unexpected operation.                 | DO NOT use the device. Call customer support for instructions on how to return the device. |

Figure 13: Error Code Display

## INSTRUMENT DEFECTS OR PRODUCT DEFICIENCIES

- For defective device concerns, or any related quality issues, please contact [info@MY01.io](mailto:info@MY01.io) or call +1 (855) 799-6901.
- A MY01 representative will deal with any quality issues related to hardware, software, or wireless functionality in a timely manner.
- Additional paper User Manuals can be made available upon request, free of charge. To request one, please email [info@MY01.io](mailto:info@MY01.io) or call +1 (855) 799-6901.

# Specifications

## DEVICE










Table 3: Device Specifications

| Device Specifications                       |                                         |
|---------------------------------------------|-----------------------------------------|
| Model                                       | MY01-0001                               |
| Pressure range                              | -99.9 to 99.9 mmHg                      |
| Power                                       | Two (2) 3 V batteries (non replaceable) |
| Battery type                                | Lithium coin cell                       |
| Display resolution                          | 0.1 mmHg                                |
| Battery life                                | 18 hours                                |
| IP Rating                                   | IP52                                    |
| Weight                                      | 85 g (+/- 5g)                           |
| Dimensions                                  | 20 cm x 6.5 cm x 3 cm (+/- 1cm)         |
| Needle gauge                                | 17-gauge                                |
| Electrical Safety                           |                                         |
| Standard                                    | EN/IEC 60601-1                          |
| Type of protection against electric shock   | Internally Powered Equipment            |
| Degree of protection against electric shock | Type BF Applied Parts                   |
| Mode of operation                           | Continuous                              |

# ENVIRONMENTAL CONDITIONS

Table 4: Environmental Conditions

|             | Operation                                                                           | Storage                                                                              | Transportation                                                                        |
|-------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Temperature |    |    |    |
| Humidity    |   |   |   |
| Pressure    |  |  |  |

## ELECTROMAGNETIC COMPATIBILITY

The MY01 device has been tested and found to comply with the requirements of EN/IEC 60601-1-2.

### WARNING

- Use of the MY01 device adjacent to or stacked with other electrical equipment should be avoided as it could result in improper operation. If such use is necessary, the MY01 device and other electrical equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MY01 device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### Note

- The MY01 device has no cables or physical accessories. Do not attempt to add any cables or accessories to the device.
- The MY01 device can communicate over a frequency band of 2.402GHz-2.4835GHz with GFSK modulation and effective radiated power of 1mW.

Table 5: Electromagnetic Emissions

| Attribute                       | Compliance/ Class Group                                              |
|---------------------------------|----------------------------------------------------------------------|
| RF Emissions Group per CISPR 11 | Group 1                                                              |
| RF Emissions Class per CISPR 11 | Class B (although the device is not for residential environment use) |

Table 6: Electromagnetic Immunity Levels

| Attribute                                                  | Compliance/ Class Group                                                                        |
|------------------------------------------------------------|------------------------------------------------------------------------------------------------|
| Electrostatic Discharge (ESD)                              | $\pm 8$ kV contact<br>$\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air<br>Per IEC 61000-4-2 |
| Radiated RF EM fields                                      | 3 V/m<br>80 MHz - 2.7 GHz<br>80 % AM at 1 kHz<br>Per IEC 61000-4-3                             |
| Proximity fields from RF wireless communications equipment | As indicated in IEC 60601-1-2                                                                  |
| Immunity to Rated Power frequency magnetic field           | 30A/m at 60Hz, 50Hz<br>Per IEC/EN 61000-4-8                                                    |

### FCC Compliance Statement

This device complies with FCC Subpart 15C rules 15.247. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

CAUTION: Any changes or modifications to this product not authorized by MY01 Inc. could void your authority to operate the product.

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This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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# Reorder information

| REF       | Product Description                            |
|-----------|------------------------------------------------|
| MY01-0001 | MY01 Continuous Compartmental Pressure Monitor |
| MY01-APP  | MY01 Mobile Application                        |

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